

Adapters Offer Safety and Durability

Offering a wide range of titanium, aluminum and carbon fiber products that are CE marked and tested to ISO 10328 standards, Freedom Innovations provides adapter and component solutions to meet the unique requirements of complex prosthetic fittings.



Adapters – General Use

	Description	Part Number	Weight Rating	Material
	Tube Clamp, 30 mm, Female Pyramid	ACC-00-12500-00	300 lbs/136 kg	Titanium
	Tube Clamp, 30 mm, Male Pyramid	ACC-00-12600-00	300 lbs/136 kg	Titanium
	Tube w/Female Pyramid, 30 mm x 100 mm	ACC-00-12700-00	300 lbs/136 kg	Titanium
	Tube w/Titanium Female Pyramid, 30 mm x 400 mm	ACC-00-12800-00	300 lbs/136 kg	Aluminum
	Pyramid, 4 Hole, Male	ACC-00-13000-00	300 lbs/136 kg	Titanium
	Pyramid, 4 Hole, Female	ACC-00-13200-00	300 lbs/136 kg	Titanium
	Male Pyramid, Rotating, 4 Hole Base	ACC-00-13500-00	300 lbs/136 kg	Titanium
	Female Pyramid, Rotating, 4 Hole Base	ACC-00-13600-00	300 lbs/136 kg	Titanium
	Base, 4 Hole, Rotating Pyramid	ACC-00-13700-00	300 lbs/136 kg	Titanium
	Base, 3 Prong, Rotating Pyramid	ACC-00-13800-00	220 lbs/99 kg	Stainless Steel
	Pyramid, Rotating, Male	ACC-00-13900-00	300 lbs/136 kg	Titanium
	Pyramid, Rotating, Female	ACC-00-14000-00	300 lbs/136 kg	Titanium





Adapters - General Use

Description	Part Number	Weight Rating	Material
 Tube, 30 mm graphite, 15", 3 mm wall, Aluminum Inserts	ACC-00-14500-00	220 lbs/99 kg	Graphite, Aluminum
 Tube, 30 mm graphite, 9", 3 mm wall, Aluminum Inserts	ACC-00-14600-00	220 lbs/99 kg	Graphite, Aluminum
 Tube, 30 mm graphite, 15", 4 mm wall, Aluminum Inserts	ACC-00-14815-00	300 lbs/136 kg	Graphite, Aluminum
 Tube, 30 mm graphite, 9", 4 mm wall, Aluminum Inserts	ACC-00-14809-00	300 lbs/136 kg	Graphite, Aluminum
 Pyramid, 4 Prong, Male	ACC-00-16000-00	300 lbs/136 kg	Titanium
 Lamination Ring, 4 Hole	ACC-00-16100-00	300 lbs/136 kg	Stainless Steel
 Adapter, Double Female Pyramid, 32 mm	ACC-00-17032-00	300 lbs/136 kg	Titanium
 Adapter, Double Female Pyramid, 45 mm	ACC-00-17045-00	300 lbs/136 kg	Titanium
 Adapter, Double Female Pyramid, 60 mm	ACC-00-17060-00	300 lbs/136 kg	Titanium
 Adapter, Double Female Pyramid, 75 mm	ACC-00-17075-00	300 lbs/136 kg	Titanium
 Alignment Jig Pyramid, Assembly with 10 mm bolt (for use with category 6-9 foot products)	ACC-00-20100-00	N/A	Titanium

Adapters – Silhouette Products

	Description	Part Number	Weight Rating	Material
	Tube Clamp, 34 mm, Female Pyramid	ACC-00-12300-00	365 lbs/165 kg	Titanium
	Tube Clamp, 34 mm, Male Pyramid	ACC-00-19600-00	365 lbs/165 kg	Titanium
	Insert to 34 mm Tube Clamp, Silhouette, Custom Fabrication	ACC-00-12400-00	365 lbs/165 kg	Nylon
	34 mm Offset Tube Clamp	ACC-00-20000-00	300 lbs/136 kg	Titanium
	Versa Clamp Adapter Kit for Pin Systems or to Increase Pylon Height (includes extended tube clamp and spacer rings)	KIT-00-123V0-00	365 lbs/165 kg	Titanium, Nylon
	Extended Tube Clamp, 34 mm, Female Pyramid	ACC-EX-12300-00	365 lbs/165 kg	Titanium
	Spacer Ring	KIT-00-12350-00	365 lbs/165 kg	Nylon
	Knee Pyramid Adapter, Male 30 mm	ACC-00-19500-00	275 lbs/124 kg	Aluminum
	Mounting Bracket Kit for Innovative Neurotronics's V-Hold™	ACC-00-20200-00	365 lbs/165 kg	Graphite

Adapters - Heavy Users

	Description	Part Number	Weight Rating	Material
	Heavy Duty Adapter Kit (includes two clamps, one 4-ear connector, one tube)	KIT-00-20200-00	500lbs/226 kg	Titanium, Graphite
	Tube Clamp 34 mm	ACC-00-20201-00	500lbs/226 kg	Titanium
	4 Ear Connector for 34 mm Sys	ACC-00-20202-00	500lbs/226 kg	Titanium
	Carbon Pylon/Tube	ACC-00-20203-00	500lbs/226 kg	Graphite