

## PARYDONN® PROSTHETIC LINER

### *The new world of user-friendly liners*

- Silky smooth—Parylene coated inner and outer surfaces feel pleasant against the skin
- No fabric covering—therefore easy to trim, no frayed edges, and easy for the amputee to keep clean and looking new
- Roll on or pull on—smooth inner surface gives amputees their choice of donning methods, a great feature for those with limited dexterity
- Just wash and go—hydrophobic nature of the coating makes the liner easy to clean and fast to dry
- Bacteria incursion blocked—virtually eliminating odors
- Reduced potential for skin complications—coating is chemical and fungal resistant
- Eliminates cold flow—proven platinum-cured silicone retains original shape even after prolonged use
- Minimized pistoning—two optimized matrix layers and varying umbrella sizes ensure a correct fit for amputees of all sizes and increase durability



PARYDONN®  
PROSTHETIC LINER

# PARYDONN<sup>®</sup> PROSTHETIC LINER

Sample Part Number  
Cushion Liner

FPL 3C-xxx

FPL 3C-220

Size in cm (22.0)

Sample Part Number  
Pin Locking Liner

FPL 3P-xxx

FPL 3P-265

Size in cm (26.5)

## Sizes: 22, 23.5, 25, 26.5, 28, 30, 32, 34 cm

Determine liner size by taking a circumferential measurement, at normal pressure, 4 cm above the distal end of the residual limb. If the exact measurement size is not available, defer to the next size down. Choosing a liner size larger than the measurement taken is not recommended.

**Wall Thickness:** 3mm

**Configurations:** Cushion or Pin Locking

**Warranty:** 6 months

**Return Policy:** 30 Days, unopened product only

