

# Plié 2.0 MPC Knee - Cosmetic Cover Kit

Cosmetically finishing the Plié 2.0 MPC Knee has never been easier. The new Plié Cosmetic Cover Kit provides a simple, snap-on solution for the lower shank portion of the prosthesis. Simply remove the Protective Cover from the knee frame, slip the knee into the pre-shaped foam and snap it into the plastic fairing that's pre-bonded to the inner core of the foam.

## Product Specifications

- Accommodates a broad range of user heights from 4'11" [151 cm, 5% female] to 6'1" [186 cm, 95% male]
- Seamless, single piece of foam simplifies the shaping process, ensures an aesthetically pleasing outcome
- The plastic fairing that is pre-bonded to the foam's inner core provides an intimate fit with the knee to prevent motion and noise
- Pre-shaped like a lower limb to reduce time spent shaping away unnecessary, excess material
- Left and right shaped foam facilitates an easy to achieve, anatomically correct shape

## Material Benefits

- Durable EVA foam, so it can take a beating
- Compared to Plastazote, EVA offers:
  - Tensile strength 2.6X greater
  - Tear strength 6.8X greater
  - Elongation at break 1.4X greater
  - More spring-like characteristics for improved durability in dynamic region near the foot/ankle transition
- Closed cell foam resists moisture and can be easily cleaned using rubbing alcohol
- Silver additives within the foam are known to provide anti-microbial and anti-fungal properties

## Part Numbers

Right: KS2-00-COVRG-00

Left: KS2-00-COVLf-00



For more information go to: [www.freedom-innovations.com/knees](http://www.freedom-innovations.com/knees)

**Global Headquarters**  
30 Fairbanks, Suite 114  
Irvine, CA 92618  
toll-free 888.818.6777  
tel 949.672.0032  
[www.freedom-innovations.com](http://www.freedom-innovations.com)

**European Head Office**  
Jaargetijdenweg 4  
7532 SX Enschede  
The Netherlands  
toll-free 00800-2806 2806  
tel +31 (0)53-20 30 300