

Microprocessor-Controlled Knee

Product	Key Features	Avg. Clearance	Avg. Weight	Weight Rating	Alignment	Warranty
Plié® 2.0 MPC Knee	Advanced stumble recovery, water resistant, user replaceable battery	9.25 in / 235 mm	1,235 g	275 lb / 125 kg high impact: 220 lb / 100 kg	Stable, 0-5 mm anterior	36 mo (12 & 24 mo service)

Carbon Fiber Prosthetic Feet

Product	Key Features	Avg. Clearance*	Avg. Weight	Weight Rating	Sizes	Warranty	LCode†
Renegade® A-T & LP-A-T	All-Terrain, Z-Shock™, EnduraCore™ Technology, lightweight design	LP: 4.7-5.6 in / 120-143 mm Tall: 6.1-7.4 in / 154-189 mm	510-515 g	365 lb / 166 kg	22-31 cm	Carbon-36 mo	L5987, L5986
Renegade® & LP	Z-Shock™ Technology, lightweight, sole plate, optional sandal-toe (22-28 cm)	LP: 4.7-5.6 in / 120-143 mm Tall: 6.1-7.4 in / 154-189 mm	510-515 g	365 lb / 166 kg	22-31 cm	Carbon-36 mo	L5987
Renegade® MX & LP-MX	Z-Shock™ Technology, with multi-axial rotation optional sandal-toe (22-28 cm)	LP: 5.6-6.5 in / 143-166 mm Tall: 6.7-8.1 in / 170-206 mm	850-945 g	285 lb / 130 kg	22-31 cm	Carbon-36 mo Rotator-24 mo	L5987, L5986
Thrive™ VS²	Load activated with Dual Keel Technology, vertical shock	5.75-7.0 in / 145-175 mm	510 g	365 lb / 166 kg	22-31 cm	36 mo	L5981, L5986 (L5999)
Runway® & HX	User adjustable heel height foot/ankle system, push button or HX design optional sandal-toe (22-28 cm)	3.6-3.9 in / 92-98 mm	485 g	255 lb / 116 kg	22-28 cm	Carbon-36 mo Ankle-18 mo	L5981, L5990
Silhouette® VS & LP-VS	Vertical shock, multi-axial, easy to cosmetically finish	LP: 6.6-7.5 in / 168-191 mm Tall: 7.0-7.9 in / 178-201 mm	300-395 g	365 lb / 166 kg	22-31 cm	36 mo	L5987, L5986
Silhouette® & LP	Lightweight, multi-axial, easy to cosmetically finish	LP: 6.6-7.5 in / 168-191 mm Tall: 7.0-7.9 in / 178-201 mm	300-395 g	365 lb / 166 kg	22-31 cm	36 mo	L5980, L5986
Promenade®	Designed for low impact users, partial split toe	4.6-5.3 in / 116-135 mm	345 g	325 lb / 147 kg	19-31 cm	36 mo	L5981
Senator®	Sophisticated simplicity provides energy return and durability, foot shell included, Senator categories 1-6	4.7-5.9 in / 120-149 mm	405 g	300 lb / 136 kg	22-30 cm	24 mo	L5976
Sierra®	Split toe design, longer keel, highly responsive	5.5-6.8 in / 140-172 mm	390 g	365 lb / 166 kg	22-31 cm	36 mo	L5980, L5986
Highlander®	Split toe design, favored by TF amputees	3.9-4.7 in / 98-121 mm	340 g	365 lb / 166 kg (HD up to 500 lb / 226 kg)	22-31 cm	36 mo	L5981, L5986
Pacifica® & LP	Low profile, split toe design, easy to cosmetically finish	LP: 1.9-2.3 in / 47-60 mm Tall: 2.3-2.5 in / 59-63 mm	320-360 g	365 lb / 166 kg	LP: 22-31 cm Tall: 24-31 cm	36 mo	L5981
LP Symes™	Lightweight, flexible, durable with direct lamination cup	1.7-2.1 in / 43-53 mm	295 g	365 lb / 166 kg	22-31 cm	36 mo	L5981
Defender®	Flexible and durable pediatric foot for users weighing up to 100 lbs	2.4-3.0 in / 61-76 mm	107 g	100 lb / 45 kg	16-21 cm	18 mo	L5981
Nitro® Running	Custom designed for runners, tall or low profile	Low 7.25 in / 184 mm Tall 9.5 in / 241 mm	276-303 g	365 lb / 166 kg	N/A	12 mo	L5980
Slalom™ Ski	Custom made; inserts directly into ski binding	5.5 in / 140 mm	1,056 g	250 lb / 114 kg	N/A	12 mo	L5981
Freestyle™ Swim	Adjust to 70° plantar flexion; push button design	4-4.6 in / 102-117 mm	542 g	255 lb / 116 kg	22-30 cm	Carbon-36 mo Ankle-12 mo	L5981, L5990
WalkTek™	iMX™ stability, optional sandal-toe (22-28 cm)	3.4-3.8 in / 85-97 mm	630 g	300 lb / 136 kg	21-30 cm	12 mo	L5972, L5986

*Clearance measurements include the foot shell and are measured from the bottom of the foot assembly to the middle of the pyramid skirt or point where the next adapter will connect. Weights are averages for 26 cm, category 4 products.
†Freedom Innovations recommends LCodes as a courtesy. The responsibility for appropriate billing lies with the attending prosthetist.

Dynamic Performance

Style and Flexibility

Everyday Performance

Unique by Design

Extreme Performance

K2 Level