

May 26, 2011

ACCESSORIES

P/N	Product Name	Weight Limit	Item	Torque Spec
AC-123	34mm Female Tube Clamp	166kg / 365lbs	Alignment set screw Pylon clamp cap screw	15Nm / 11ft-lb 15Nm / 11ft-lb
AC-123L	Extended Tube Clamp, 34mm Female	166kg / 365lbs	Alignment set screw Pylon clamp cap screw	15Nm / 11ft-lb 15Nm / 11ft-lb
AC-125	30mm Female Tube Clamp	136kg / 300lbs	Alignment set screw Pylon clamp cap screw	15Nm / 11ft-lb 5Nm / 4ft-lb
AC-126	30mm Male Tube Clamp	100kg / 220lbs	Alignment set screw	7Nm / 5ft-lb
AC-127	30 x 100mm Female Tube	136kg / 300lbs	Alignment set screw	15Nm / 11ft-lb
AC-128	30 x 400mm Female Tube	100kg / 220lbs	Alignment set screw	15Nm / 11ft-lb
AC-130	4-Hole Male Pyramid	136kg / 300lbs	Attachment screw	15Nm / 11ft-lb
AC-132	4-Hole Female Pyramid	136kg / 300lbs	Alignment set screw Mounting screw	15Nm / 11ft-lb 15Nm / 11ft-lb
AC-135	4-Hole Rotating Male Pyramid	136kg / 300lbs	Alignment & mounting screw Rotation clamp screw	15Nm / 11ft-lb 5Nm / 4ft-lb
AC-136	4-Hole Rotating Female Pyramid	136kg / 300lbs	Alignment & mounting screw Rotation clamp screw	15Nm / 11ft-lb 5Nm / 4ft-lb
AC-137	4-Hole Rotating Pyramid Base	136kg / 300lbs	Mounting screw Rotation clamp screw	15Nm / 11ft-lb 5Nm / 4ft-lb
AC-138	3-Prong Rotating Pyramid	100kg / 220lbs	Mounting screw Rotation clamp screw	15Nm / 11ft-lb 5Nm / 4ft-lb
AC-139	Rotating Male Pyramid	136kg / 300lbs	NA	NA
AC-140	Rotating Female Pyramid	136kg / 300lbs	Rotation clamp screw	5Nm / 4ft-lb
AC-145	30mm x 3mm x 15in Composite Tube	100kg / 220lbs	NA	NA

AC-146	30mm x 3mm x 9in Composite Tube	100kg / 220lbs	NA	NA
AC-147	Anodized Inserts for 3mm Composite Tube	100kg / 220lbs	NA	NA
AC-148	30mm x 4mm x 15in Composite Tube	136kg / 300lbs	NA	NA
AC-149	30mm x 4mm x 9in Composite Tube	136kg / 300lbs	NA	NA
AC-150	Anodized Inserts for 4mm Composite Tube	136kg / 300lbs	NA	NA
AC-160	4-Prong Male Pyramid	136kg / 300lbs	Mounting screw	15Nm / 11ft-lb
AC-161	4-Hole Lamination Base	136kg / 300lbs	Mounting screw	5Nm / 4ft-lb
AC-170	32mm Dual Female Adaptor	136kg / 300lbs	Alignment set screw	15Nm / 11ft-lb
AC-171	45mm Dual Female Adaptor	136kg / 300lbs	Alignment set screw	15Nm / 11ft-lb
AC-172	60mm Dual Female Adaptor	136kg / 300lbs	Alignment set screw	15Nm / 11ft-lb
AC-173	75mm Dual Female Adaptor	136kg / 300lbs	Alignment set screw	15Nm / 11ft-lb
AC-195	Knee Pyramid Adaptor	136kg / 300lbs	NA	NA
AC-196	34mm Male Tube Clamp	166kg / 365lbs	Pylon clamp cap screw	15Nm / 11ft-lb
AC-200	Offset Tube Clamp, 34mm Female	166kg / 365lbs	Alignment set screw Pylon clamp cap screw	15Nm / 11ft-lb 15Nm / 11ft-lb
AC-202-CLA	34mm Female Tube Clamp	226kg / 500lbs	Alignment set screw Pylon clamp cap screw	15Nm / 11ft-lb 18Nm / 13ft-lb
AC-202-EAR	4-Prong Male Pyramid	226kg / 500lbs	Alignment/mounting screw	15Nm / 11ft-lb
AC-202-PYL	34mm Composite Pylon	226kg / 500lbs	NA	NA

REVOLUTION SERIES

P/N	Product Name	Weight Limit	Item	Torque Spec
RS1	Renegade 19-21 cm (Regular toe only)	45kg / 100lbs	M10 Connector Bolt (Cat 1-6)	40-45 ft-lbs & Loctite
	Renegade Mold 22-31 cm (Regular & Sandal)	166kg / 365lbs	M10 Connector Bolt (Cat 1-6) M12 Connector Bolt (Cat 7-9)	40-45 ft-lbs & Loctite 45-50 ft-lbs & Loctite
RS1-M1	Renegade MX	Cat 1-7 130kg / 285lbs	NA	NA
RS2-P1	Runway Push Button	Cat 1-6 166kg / 255lbs	Slide Bolts (Cat 1-6)	17 ft-lbs & Loctite
RS2-K1	Runway HX	Cat 1-6 166kg / 255lbs	Slide Bolts (Cat 1-6)	17 ft-lbs & Loctite
RS3	Renegade LP (International) 22-25 cm	166kg / 365lbs	M10 Connector Bolt (Cat 1-4) M12 Connector Bolt (Cat 5-9)	40-45 ft-lbs & Loctite 45-50 ft-lbs & Loctite
	Renegade LP (International) 26-28 cm	166kg / 365lbs	M10 Connector Bolt (Cat 1-6) M12 Connector Bolt (Cat 7-9)	40-45 ft-lbs & Loctite 45-50 ft-lbs & Loctite
	Renegade LP (International) 29-31 cm	166kg / 365lbs	M12 Connector Bolt (Cat 1-9)	45-50 ft-lbs & Loctite
RS4	Silhouette - 3/8" Heel	166kg / 365lbs	NA	NA
RS5	Silhouette LP - 3/8" Heel	166kg / 365lbs	NA	NA
RS6	Reneage LP (US) 22-25 cm	166kg / 365lbs	M10 Connector Bolt (Cat 1-6) M12 Connector Bolt (Cat 7-9)	40-45 ft-lbs & Loctite 45-50 ft-lbs & Loctite
	Reneage LP (US) 26-28 cm	166kg / 365lbs	M10 Connector Bolt (Cat 1-4) M12 Connector Bolt (Cat 5-9)	40-45 ft-lbs & Loctite 45-50 ft-lbs & Loctite
	Reneage LP (US) 29-31 cm	166kg / 365lbs	M12 Connector Bolt (Cat 1-9)	45-50 ft-lbs & Loctite
RS6-M1	Renegade LP MX	Cat 1-7 130kg / 285lbs	NA	NA
RS8	Silhouette - 3/4" Heel	166kg / 365lbs	NA	NA
RS9	Silhouette LP - 3/4" Heel	166kg / 365lbs	NA	NA
R14	Silhouette VS- 3/8" Heel	166kg / 365lbs	NA	NA
R15	Silhouette VS LP - 3/8" Heel	166kg / 365lbs	NA	NA
R18	Silhouette VS - 3/4" Heel	166kg / 365lbs	NA	NA
R19	Silhouette VS LP - 3/4" Heel	166kg / 365lbs	NA	NA

FREEDOM SERIES

P/N	Product Name	Weight Limit	Item	Torque Spec
FS1	Sierra 19-21 cm	45kg / 100 lbs	M10 Connector Bolt (Cat 1-6)	40-45 ft-lbs & Loctite
	Sierra 22-31 cm	166kg / 365lbs	M10 Connector Bolt (Cat 1-6) M12 Connector Bolt (Cat 7-9)	40-45 ft-lbs & Loctite 45-50 ft-lbs & Loctite
FS2	Pacifica 22-31 cm	166kg / 365lbs	NA	NA
FS3	Defender 16-18 cm	Cat 1P-4A 45kg / 100lbs	M8 Connector Bolt (Cat 1P-1A) M10 Connector Bolt (Cat 2A-4A)	16-18 ft-lbs & Loctite 40-45 ft-lbs & Loctite
	Defender 19-21 cm	Cat 1P-4A 45kg / 100lbs	M8 Connector Bolt (Cat 1P-1A) M10 Connector Bolt (Cat 2A-4A)	16-18 ft-lbs & Loctite 40-45 ft-lbs & Loctite
FS3	Highlander 22-25 cm	166kg / 365lbs	M10 Connector Bolt (Cat 1-6) M12 Connector Bolt (Cat 7-9)	40-45 ft-lbs & Loctite 45-50 ft-lbs & Loctite
	Highlander 26-28 cm	166kg / 365lbs	M10 Connector Bolt (Cat 1-4) M12 Connector Bolt (Cat 5-9)	40-45 ft-lbs & Loctite 45-50 ft-lbs & Loctite
	Highlander 29-31 cm	166kg / 365lbs	M12 Connector Bolt (Cat 1-9)	45-50 ft-lbs & Loctite
FS3	Highlander MAX	227kg / 500lbs	M12 Connector Bolt (Cat 10-13)	45-50 ft-lbs & Loctite
FS4	Pacifica LP 22-25 cm	166kg / 365lbs	M8 Pyramid Bolts (Cat 1-6) M10 Pyramid Bolts (Cat 7-9)	15 ft-lbs & Loctite 20 ft-lbs & Loctite
	Pacifica LP 26-28 cm	166kg / 365lbs	M8 Pyramid Bolts (Cat 1-4) M10 Pyramid Bolts (Cat 5-9)	15 ft-lbs & Loctite 20 ft-lbs & Loctite
	Pacifica LP 29-31 cm	166kg / 365lbs	M10 Pyramid Bolts (Cat 1-9)	20 ft-lbs & Loctite
	Thrive	166kg / 365lbs	M12 Connector Bolt (Cat 1-9)	45-50 ft-lbs & Loctite
FS6	Nitro	166kg / 365lbs	M10 Connector Bolt (Cat 1-4)	40-45 ft-lbs & Loctite
			M12 Connector Bolt (Cat 5-9)	45-50 ft-lbs & Loctite
FS9	Nitro LP	166kg / 365lbs	M10 Connector Bolt (Cat 1-4)	40-45 ft-lbs & Loctite
			M12 Connector Bolt (Cat 5-9)	45-50 ft-lbs & Loctite

SPECIALTY SERIES

P/N	Product Name	Weight Limit	Item	Torque Spec
VS1	Senator	Cat 1-6* 136kg / 300lbs	M10 Connector Bolt (Cat 1-6)*	40-45 ft-lbs & Loctite
VS2	Promenade	Cat 1-5* 148kg / 325lbs	M10 Connector Bolt (Cat 1-5)*	40-45 ft-lbs & Loctite
VS3	WalkTek	Cat 1-4* 136Kg/300lbs	NA	NA
LP2	LP Symes	166kg / 365lbs	M10 Cup Bolts (Cat)	35ft-lbs
LP2-W2	Swim	Cat 1-6 116kg / 255lbs	M10 Bolts	7 ft-lbs & Loctite
SKI	Ski Foot	113kg / 250lbs	M12 Cup Bolts M8 Toe Bolts	45-50 ft-lbs & Loctite 25-30 ft-lbs & Loctite
* Special category scale applies				









